mē[̄] my eļōs

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Hello and congratulations on obtaining your new mē facial kit! What is this kit?

This kit allows you to safely remove hair from your face with mē. A booklet containing important information and protective eyewear has been included in order for you to safely use the mē device on your face.

Removing unwanted hair from the face is very different from removing unwanted hair from other parts of the body. Please be sure to read this booklet thoroughly before starting to remove facial hair with mē.

Be sure to read the mē user guide, which includes detailed instructions and precautions, before starting treatment anywhere with the mē device!

Why is hair removal from your face not like other body areas?

- The hair on the face differs from other body areas in its physiology. Its growing cycles are different, and its characteristics, such as thickness (vellus hair), color and depth of follicle are different from body hair.
- The face is more sensitive to hair removal due to sensitivity of nerves in some of its areas.
- The skin on the face and neck is more likely to be exposed daily to the sun's UV radiation and therefore may be more sensitive than other parts of your body.
- Treating near your eyes can cause temporary discomfort, and therefore it is required to use protective eyewear and treat only from your cheeks and down.

When you should not use mē on the face

- If your skin type is 4-6 (approved for skin types 1-3 only).

- If your facial skin is sunburned or tanned, wait until the sunburn or tan fades prior to performing a test pulse on a hidden area. Evaluate the area after 3 days to ensure no side effects have occurred and only then commence treatment.
- If you have undergone any professional facial treatment (such as chemical peel, laser/light/radio frequency treatment or injected dermal fillers) within the past 3 months.

- If you have undergone any facial surgery (such as face lift, facial implant) within the last three months or before complete healing, please consult with your doctor.
- If you are taking any medications or have a history of any skin disorders, consult with your doctor before using mē on the face or any other body area.
- Please refer to mē User Guide for general list of exclusion criteria.

Precautions

- Treat below the Cheeks line only.

- DO NOT use the epilator or shaver accessories when using mē on the face.
- DO NOT sun tan or use sun beds for at least 2 weeks before a treatment.
- Avoid starting a treatment on your face when you know you will be exposing your face to sunlight in the following 24 hours. If you are going into the sun following a treatment take care to shade the treated area from direct sunlight using a hat or parasol.
- After treatment avoid exposing the treated area to the sun for at least 24 hours.

- Apply sunscreen of at least SPF 30 on the treated area for 2 weeks following the treatment.
- DO NOT treat the eye socket or around the eye area.
- DO NOT use for eyebrow shaping! It is important to avoid accidentally treating the eyebrows as any hair loss is likely to be permanent.
- DO NOT treat over moles or birth marks.
Protect small areas of sensitive skin such as moles, birthmarks etc by using a white eyeliner.
For larger areas of sensitive skin, use a piece of card or thick paper to cover the affected area and prevent exposure to the elōs treatment.
- Protect lips and the lip margin using a white eyeliner.
- DO NOT treat tattoos or permanent makeup.
- DO NOT treat the scalp.
- DO NOT treat inside the nostrils or ears as these areas are delicate.
- Avoid treating areas that contain only vellus hair

(very thin, fine hair, also called “peach fuzz”), as this may result in undesired hair re-growth.

- In the very unlikely event that you see pigmented changes on your face please stop the treatment and consult with your doctor.
- Do not use perfumed products or make up, for up to 24 hours after a treatment as this may cause slight irritation.
- If unusual pain occurs during treatment, you should stop treatment.
- Temporary skin redness may occur up to 24 hours after treatment. If the symptoms persist more than 24 hours, and/or you feel significant discomfort and/or encounter any severe side effects, discontinue treatment and consult with your doctor.

How to use

1. The treatment should be performed in a broad ambient light environment.
2. Adjust the eyewear's elastic band so they fit comfortably over your eyes.
3. Conceal the delicate areas such as lip edges and facial moles with a white eyeliner.
4. Turn on your mē device.
5. Select the LOW elōs level on the base unit and press the elōs button on your applicator.
6. Place the protective eyewear so that it covers your eyes.

7. Apply a single evaluation pulse on a hidden area on your face or neck and observe treated area (If you observe any unusual reaction, stop the treatment immediatly).
8. Hold the upper part of your mē applicator without touching the two grip sensors and gently position it over the area you want to treat.
9. To allow elōs flashes on the area you want treat, touch momentarily both grip sensors with the thumb and forefinger of your other hand and release after 1-2 pulses were emitted.

grip sensors

10. Without holding the grip sensors, Reposition the applicator on the next area to be treated and repeat step 8.

* Always start with LOW elōs level. If you feel comfortable, you may increase the elōs level for enhanced efficacy.

Tips

- Please take additional care when treating hair above the upper lip because the central area below the nostrils is particularly sensitive.
- Hair growth under the chin can be denser which can lead to light reddening in the skin following a treatment. This redness is temporary and should completely disappear after 24 hours.
- After the treatment we recommend using an oil free moisturizer for the face in order to feel smoothed.

- For men it is not advisable to treat the beard as it might create uneven results, which may produce changes in facial characteristics that may not be desirable.

For more information regarding facial treatment with mē device, please visit www.memyelos.com

The Warranty granted by Manufacturer and Distributor to the mē product, included within the mē packaging, shall be the sole and exclusive warranty applicable to this Facial Kit. The Disclaimer of Warranties in the Warranty shall also apply to this Facial Kit.

Distributed by HOMEDICS GROUP LTD
PO Box 460, Tonbridge, Kent TN9 9EW, UK
Tel: 0203 302 0623 Fax: +44 (0) 131 202 0697
E-mail: meservice@homedics.co.uk Website: www.homedics.co.uk

Manufactured by Syneron Medical Ltd
P.O. Box 550 Industrial Zone, Yokneam Illit 20692, Israel
Tel: +972 73 2442200 Website: www.memyelos.com

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© 2011. All rights reserved. Syneron, the Syneron Logo, mē, my elōs, and elōs are trademarks of Syneron Medical Ltd. and may be registered in certain jurisdictions. elōs (electro-optical synergy) is a proprietary technology of Syneron Medical. Patent pending.
Hergestellt von Syneron Medical Ltd